



ACUTE CORONARY SYNDROMES

Clinical Utility of Cardiac Troponin



PIONEERING DIAGNOSTICS



PREFACE

Coronary heart disease (CHD) remains one of the major causes of death worldwide. Presentation with CHD, typically as chest pain, is one of the more common conditions seen in the Emergency Department (ED).

- Although chest pain is a common complaint, patients who have myocardial infarction (MI) are a relatively small percentage of all those who present with chest pain. Early recognition of those without MI is therefore as important as early identification of those who have MI.
- New developments in measurement methods mean that cardiac troponin (cTn) can be measured reliably and repeatedly at much lower levels than was previously possible. These methods are called high-sensitivity cardiac troponin (hs-cTn) assays.
- The use of hs-cTn measurement means it is possible to exclude myocardial infarction by measurement on admission to hospital and as little as 1-2 hours later.
- Rapid diagnostic pathways mean that both long waits in the ED while test results are awaited can be avoided and unnecessary hospital stays shortened.

This booklet provides clinicians and laboratory managers with concise, up-to-date information on the pathophysiology and diagnosis of myocardial infarction with special emphasis on the role of hs-cTn measurements for risk stratification and management of patients. Aspects of the laboratory and clinical factors affecting the routine measurement of cTn using hs-cTn assays are also discussed.

OUR SPECIAL THANKS GO TO

Professor Paul Collinson

MA, MB, BChir, FRCPath, MD, FACB, FRCP
Consultant Chemical Pathologist,
St George's University Hospitals NHS Foundation Trust;
Professor of Cardiovascular Biomarkers,
St. George's University of London, UK

Professor Yann-Erick Claessens

MD, PhD, MBA
Head of Emergency Department,
Princess Grace Hospital, Monaco;
Professor of Emergency Medicine,
Université Paris Descartes, Paris, France



Professor Paul Collinson



Professor Yann-Erick Claessens

ACUTE CORONARY SYNDROMES: LEADING CAUSE OF DEATH

High mortality and morbidity

- Globally, the total number of deaths due to **ischemic heart disease** (IHD) has risen steadily to reach 9.14 million in 2019.¹ This makes IHD the leading cause of death in the world. Over the same time period the global age-standardized IHD death rate declined, indicating that, on average, global increase in IHD deaths has been due to population growth and aging.
- In 2019, there were an estimated 197 million prevalent cases of IHD worldwide.¹ In the 54 member countries of the European Society of Cardiology (ESC) alone there were 34.9 million people living with IHD in 2017.² The societal impact of IHD is huge and the global burden of IHD, in terms of disability-adjusted life years (DALY's), has increased since 1990 with IHD the leading cause of DALY's in people 50 years and older.³ The DALY, a measure of the burden of a disease in populations, combines life expectancy (years of life lost) with adjusted quality of life (years lived with the disease).

Chest pain and acute coronary syndromes

- **Chest pain** is the chief complaint in patients with **cardiac ischemia**.
- **Chest pain** is among the top-10 reasons to visit an **emergency department** (ED), accounting for about 5% of all ED visits.⁴
- Among patients presenting with acute chest pain to the ED, disease prevalence can be expected as follows:⁵
 - **Acute coronary syndromes (ACS) in 35%**
 - ST-segment elevation myocardial infarction (**STEMI**) in 5-10%
 - Non-ST-segment elevation myocardial infarction (**NSTEMI**) in 15-20%
 - Unstable angina (**UA**) in 10%
 - Other cardiac conditions in 15%
 - Non-cardiac diseases in 50%
- **Acute myocardial infarction** (AMI: STEMI plus NSTEMI)* accounts for about 2/3rd and UA for about 1/3rd of ACS. There is a global trend of a decreasing incidence in STEMI with a concomitant increase in NSTEMI.⁶ STEMI is relatively more common in younger than in older people, and is more common in men than in women.⁶

* The **Universal Definition of Myocardial Infarction** (currently in its 4th update) combines the ECG, clinical features and cardiac troponin (cTn) values above the 99th percentile upper reference limit (URL) with cTn changes (cTn kinetics, delta values) as an integral part of the definition.⁷ The shift to cTn testing has had an impact on the epidemiology of AMI by reclassifying patients previously diagnosed as unstable angina to NSTEMI.

CONTENTS

ACUTE CORONARY SYNDROMES

- 1 Definition, classification and pathophysiology p. 4
- 2 Diagnosis p. 10

CARDIAC TROPONIN

- 1 Biochemistry and pathophysiology p. 14
- 2 Assays p. 15
- 3 Role of cardiac troponin in diagnosis of ACS p. 16
- 4 Guide to interpretation p. 22

FREQUENTLY ASKED QUESTIONS p. 24

REFERENCES p. 32

LIST OF ABBREVIATIONS & ACRONYMS p. 36

➤ For easy reading and reference, look for the blue boxes highlighting the key points in each chapter.

ACUTE CORONARY SYNDROMES

1 Definition, classification and pathophysiology

The term **acute coronary syndromes (ACS)** refers to a constellation of clinical signs and symptoms due to acute **myocardial ischemia** that results from inadequate delivery of oxygen to the heart muscle. In myocardial ischemia there is an imbalance between myocardial oxygen supply and/or demand. This usually occurs due to reduced blood flow because of partial or complete blockage of the coronary arteries.


ACS may range from a potentially reversible phase (**unstable angina, UA**) to irreversible cell death (**myocardial infarction, MI**). The diagnosis of ACS requires the integration of clinical signs and symptoms, ECG and serial cardiac troponin (cTn) measurements (**Figure 1**).

- The **ECG** identifies approximately one quarter of ACS patients who present with persistent (>20 min) ST-segment elevation. These patients are diagnosed with **STEMI** and require immediate reperfusion.⁶
- The other three quarters of ACS patients are without persistent ST-segment elevations (NSTEMI-ACS; non-ST-elevation ACS). In these patients **measurement of cTn** is required to distinguish NSTEMI from UA (where there is myocardial ischemia in the absence of myocardial necrosis).⁸

The distinction between STEMI, NSTEMI and UA is clinically important and drives the decision for timing, type and intensity of therapeutic intervention.^{6, 8}

Figure 1. Diagnostic algorithm and triage in ACS.⁸

Adapted from Collet JP, et al. *Eur Heart J*. 2021;42:1289-1367



*Noncardiac: thoracic diseases such as pneumonia or pneumothorax.

**Other cardiac: including amongst others, myocarditis, Takotsubo syndrome, congestive heart failure.

Myocardial infarction and myocardial injury

The 4th Universal Definition of Myocardial Infarction introduces a differentiation of myocardial infarction from myocardial injury.⁷

The term **myocardial injury** should be used when there is evidence of elevated cTn values with at least one value above the 99th percentile URL. Myocardial injury is considered acute if there is a rise and/or fall of cTn values. Myocardial injury, often acute, occurs in a wide range of conditions and is three times more common in the ED than a cTn elevation due to myocardial infarction (see below).

The term **acute myocardial infarction** (types 1, 2 and 3 MI) should be used when there is acute myocardial injury with clinical evidence of acute myocardial ischemia and with detection of a rise and/or fall of cTn values with at least one value above the 99th percentile URL and at least one of the following⁷:

- Symptoms of myocardial ischemia.
- New ischemic ECG changes.
- Development of pathological Q waves.
- Imaging evidence of new loss of viable myocardium or new regional wall motion abnormality in a pattern consistent with an ischemic etiology.
- Identification of a coronary thrombus by angiography or autopsy (not for types 2 or 3 MIs).

Depending on the underlying pathophysiological mechanism, 5 types of MI are distinguished including acute MI (types 1,2 and 3 MI) and procedure-related MI (types 4 and 5 MI) (**Table 1**).

Clinical criteria for MI

The clinical definition of MI denotes the presence of acute myocardial injury detected by abnormal cardiac biomarkers in the setting of evidence of acute myocardial ischemia.

Table 1. Etiological classification of myocardial infarction.⁷

Adapted from Thygesen K, et al. *Circulation*. 2018;138:e618-e651

Type 1 MI	MI caused by atherothrombotic coronary artery disease (CAD) and usually precipitated by atherosclerotic plaque disruption (rupture or erosion).
Type 2 MI	MI due to ischemic myocardial injury and evidence of an imbalance between myocardial oxygen supply and/or demand unrelated to coronary thrombosis
Type 3 MI	MI resulting in sudden cardiac death with symptoms suggestive of myocardial ischemia but without definitive cardiac biomarker evidence of MI (patients die before blood samples for cTn testing can be obtained or before increases in cTn can be identified).
Type 4 MI	MI associated with percutaneous coronary intervention (PCI; type 4a), stent thrombosis (type 4b) or restenosis (type 4c).
Type 5 MI	MI related to coronary artery bypass grafting (CABG).

The majority of MI patients present with **type 1 MI**, but **type 2 MI** is also a common subtype that may account for up to a third of all MI cases.⁹ Patients with type 2 MI are more frequently female and are older with more comorbidities; consequently they have an almost 2-fold poorer long-term survival than type 1 MI, with excess deaths due to non-cardiovascular causes.¹⁰


The distinction between type 1 and type 2 MI and acute myocardial injury is clinically important because treatment is different.

Type 1 MI is caused by **plaque rupture** or **erosion** with **atherothrombosis** (**Figure 2**) and is treated by revascularization procedures and antithrombotic therapy.^{6,8} Plaque rupture is the most common cause both in STEMI and NSTEMI, but plaque erosion is more frequent in NSTEMI patients.¹¹

- **Plaque rupture** is the most common mechanism; it is associated with rupture of the thin fibrous cap of lipid-rich atherosclerotic plaques and exposure of the necrotic core to the vessel lumen leading to thrombus formation.¹¹
- **Plaque erosion** involves thrombus formation in an area of endothelial denudation adjacent to an atherosclerotic plaque without disruption of the fibrous cap.¹¹

Figure 2. Myocardial infarction type 1.⁷


Adapted from Thygesen K, et al. *Circulation*. 2018;138:e618-e651



Type 2 MI is the result of a **mismatch** between **myocardial oxygen supply and demand** in the **absence of atherothrombosis** (Figure 3) and consequently its treatment is based on the diagnosis of the underlying disorder.⁹ Underlying causes of type 2 MI are listed in **Table 2**.

Figure 3. Myocardial infarction type 2.⁷

Adapted from Thygesen K, et al. *Circulation*. 2018;138:e618-e651



MINOCA

It is increasingly recognized that there is a group of MI patients with no angiographic obstructive CAD, and the term *myocardial infarction with non-obstructive coronary arteries (MINOCA)* has been coined for this entity.⁷ Like the diagnosis of MI, the diagnosis of MINOCA indicates that there is an ischemic mechanism responsible for myocyte injury. Among patients diagnosed with MI, the prevalence of MINOCA is estimated to be 6- 8%; more common in women and in patients presenting with NSTEMI.⁷

Acute myocardial injury

The current update of the Universal Definition of Myocardial Infarction recognizes that myocardial injury may occur due to other pathologies than ischemic heart disease.⁷ Acute myocardial injury is suspected when there is a cTn elevation without symptoms of myocardial ischemia, or obvious CAD, or another clinical cardiologic diagnosis. It usually accompanies acute illness and requires treatment of the underlying medical condition. When it occurs, however, it indicates an adverse prognosis. It has been found that in the ED population, cTn elevation may be more common due to myocardial injury than myocardial infarction.¹² Conditions associated with non-ischemic myocardial injury are listed in **Table 2**.

Table 2. Causes of myocardial injury with elevation of cTn.⁷

Adapted from Thygesen K, et al. *Circulation*. 2018;138:e618-e651

A. Myocardial injury related to acute myocardial ischemia

Atherosclerotic plaque disruption with thrombosis (type 1 MI)

Oxygen supply/demand imbalance (type 2 MI)

Reduced myocardial perfusion, e.g.

- Coronary artery spasm, microvascular dysfunction
- Coronary embolism
- Coronary artery dissection
- Sustained bradyarrhythmia
- Hypotension or shock
- Respiratory failure
- Severe anemia

Increased myocardial oxygen demand, e.g.

- Sustained tachyarrhythmia
- Severe hypertension with or without left ventricular hypertrophy

B. Other (non-ischemic) causes of myocardial injury

Cardiac conditions, e.g

- Heart failure
- Myocarditis
- Cardiomyopathy (any type)
- Takotsubo syndrome
- Coronary revascularization procedure
- Cardiac procedure other than revascularization
- Catheter ablation
- Defibrillator shocks
- Cardiac contusion

Systemic conditions, e.g.

- Sepsis, infectious disease (including COVID-19)
- Chronic kidney disease
- Stroke, subarachnoid hemorrhage
- Pulmonary embolism, pulmonary hypertension
- Infiltrative diseases, e.g. amyloidosis, sarcoidosis
- Chemotherapeutic agents
- Critically ill patients
- Strenuous exercise

2 Diagnosis

The diagnosis of ACS (Figure 1) is based on the integration of:

- Clinical presentation (history, symptoms, signs, location and severity of pain)
- 12-lead ECG
- cTn measurement

Clinical presentation

Demographic and clinical features such as older age, male gender, family history of CAD, diabetes, hyperlipidemia, hypertension, renal insufficiency, previous manifestation of CAD as well as peripheral or carotid artery disease increase the likelihood of ACS.^{5,8} Conditions that may exacerbate or precipitate ACS include anemia, infection, inflammation, fever, and metabolic or endocrine (thyroid) disorders.⁸

The typical **chief complaint** in ACS is **angina** (Table 3), defined as central chest pain or discomfort that occurs due to inadequate delivery of oxygen to the heart muscle. It has been likened to having a heavy weight pressing down inside the chest.

The pain may radiate to the neck, jaw or arm (usually the left). Typical stable angina (often called effort angina or angina pectoralis) is provoked by exertion or emotional stress and is relieved by rest or sublingual nitrates. This does not occur with unstable angina or myocardial infarction.

Often, the discomfort is diffuse (not localized, not positional) and may be accompanied by sweating, shortness of breath (dyspnea), nausea or fainting (syncope) (Table 3). Chest pain is not specific for ACS and may also occur in many other cardiac and non-cardiac conditions (Table 4). Like ACS, some of these conditions are also life-threatening, such as acute aortic dissection (AAD) and pulmonary embolism (PE).

Unfortunately, in patients with suspected ACS, elements of chest pain history alone are highly unlikely to confirm or exclude the diagnosis of ACS due to their limited sensitivity and specificity.^{13,14}

Table 3. Clinical presentation in suspected ACS.⁸

Adapted from Collet JP, et al. *Eur Heart J*. 2021;42:1289-1367

Typical	<ul style="list-style-type: none"> • Chest pain: retro-sternal pressure or heaviness ('angina') radiating to left arm, neck or jaw, which may be intermittent or persistent. • May be accompanied by sweating, nausea, abdominal pain, dyspnea or syncope.
Atypical	<ul style="list-style-type: none"> • Epigastric pain, recent onset indigestion, stabbing chest pain, chest pain with pleuritic features or increasing dyspnea. • Often observed in younger (25-40 years) and older (>75 years) patients, in women, and in patients with diabetes, chronic renal failure or dementia.

Table 4. Differential diagnosis: causes of chest pain other than ACS.⁸

Adapted from Collet JP, et al. *Eur Heart J*. 2021;42:1289-1367

Cardiac	Gastro-intestinal
Myopericarditis Cardiomyopathies Tachyarrhythmias Acute heart failure Hypertensive emergencies Aortic valve stenosis Takotsubo syndrome Coronary spasm Cardiac trauma	Esophagitis, reflux or spasm Peptic ulcer, gastritis Pancreatitis Cholecystitis
Pulmonary	Orthopedic
Pulmonary embolism (Tension)-Pneumothorax Bronchitis, pneumonia Pleuritis	Musculoskeletal disorders Chest trauma Muscle injury/inflammation Costochondritis Cervical spine pathologies
Vascular	Other
Aortic dissection Symptomatic aortic aneurysm Stroke	Anxiety disorders Herpes zoster Anemia

Bold: common and/or important differential diagnoses.

ECG

The 12-lead ECG is the first-line diagnostic tool in the assessment of patients with suspected ACS (**Figure 1**). It should be recorded within 10 min of the patient's arrival in the ED or at first prehospital medical contact.^{6,8}

The ECG may be normal in more than one-third of suspected ACS patients. Characteristic **ECG manifestations of acute myocardial ischemia** include **ST-depression, ST-elevation** and **T-wave changes**.⁷

The ECG may also provide information on the anatomic location, extent and severity of the coronary lesion as well as the presence of complications of acute MI with a poor prognosis (e.g. arrhythmias and conduction abnormalities).

The ECG diagnosis of AMI can be complicated in the presence of left bundle branch block (LBBB) or ventricular pacing.⁶ It is also important to be aware of conditions that mimic myocardial ischemia or MI on ECG.⁷ Examples of such false-positive ECG abnormalities include cardiomyopathy, cardiac amyloidosis, LBBB, ventricular hypertrophy, myocarditis, or hyperkalemia.

If the standard leads are inconclusive and the patient has signs or symptoms suggestive of ongoing myocardial ischemia, additional leads should be recorded.⁸

Cardiac troponin


Cardiac troponin (I or T) is the preferred biomarker of cardiac necrosis.⁷

Latest 2020 ESC Guidelines for the management of ACS in patients presenting without persistent ST-segment elevation state that, for diagnostic purposes, it is not recommended to routinely measure additional biomarkers such as CK, CK-MB, h-FABP, or copeptin, in addition to hs-cTn.⁸

cTn measurement needs to be integrated with clinical assessment and ECG in the diagnosis of patients with suspected ACS (**Figure 1**). See next chapter for further details on its use in diagnostic algorithms.

Figure 4. ECG tracings in ACS.

Adapted from www.thrombosisadviser.com/acute-coronary-syndrome



Persistent ST-segment elevation is the hallmark finding in STEMI. ST-depression and T-wave inversion are notable signs of cardiac ischemia in NSTEMI.

CARDIAC TROPONIN

1 Biochemistry and pathophysiology


The troponin complex consists of 3 subunits (I, T and C) and is essential for the regulation of skeletal and cardiac muscle contraction.¹⁵ In contrast to troponin C, cardiac-specific isoforms of troponin T and I exist.

The release of cTnI and cTnT from cardiomyocytes is highly specific for myocardial injury (Figure 5). However, any type of injury, not just ischemic injury, can result in the release of cTn into the blood circulation.⁷ Consequently, cTn can also be elevated in other conditions than MI (Table 2).

Figure 5. Troponin is a marker of cardiac damage.^{15, 17}

Adapted from Parmacek MS, Solaro RJ. *Prog Cardiovasc Dis.* 2004;47:159-176; Bates KJ, et al. *Clin Chem.* 2010;56:952-958

CARDIAC MYOCYTE



The troponin complex is the regulatory complex of the myofibrillar thin filament (consisting of actin and tropomyosin) that plays a critical role in regulating Ca²⁺-mediated excitation-contraction coupling in the heart. There is one troponin complex bound to every seventh actin monomer. The troponin complex consists of 3 subunits: troponin C (TnC, 18 kDa), the Ca²⁺-binding subunit; troponin I (TnI, 23 kDa), the Inhibitory subunit that prevents contraction in the absence of Ca²⁺-binding to TnC; and troponin T (TnT, 35 kDa), the subunit that attaches troponin to tropomyosin and to the myofibrillar thin filament. The cytosol of the cardiac myocyte contains free unbound troponins which are released into the circulation upon damage. After AMI, cTnI is present in serum as ternary T-I-C and binary I-C complexes, whereas cTnT is present as a combination of T-I-C and free cTnT.


The classical viewpoint is that cTn release from cardiomyocytes is caused by irreversible injury resulting in cardiac necrosis. However, the release of cTn does not always require myocardial cell death and there are clinical scenarios with cTn increase despite a low clinical likelihood of myocardial necrosis.¹⁶ Examples include extraordinary endurance exercise (e.g. marathon), after reversible ischemia in an exercise or dobutamine stress test, after rapid atrial pacing, in severe chronic stable heart failure and in severe chronic stable artery disease.

2 Assays

Specific antibodies have been raised against the cardiac-specific cTn isoforms, the basis of widely available quantitative cTnI and cTnT assays.¹⁸ Advances in technology have led to an improvement of the ability of cTn assays to detect and quantify cardiomyocyte injury. Because of their increased analytical power, allowing accurate measurement below the 99th percentile URL as well as accurate measurement of small absolute changes, **high-sensitivity cTn (hs-cTn) assays enable earlier detection of AMI** with shortening of the time window for serial testing (Figure 6).¹⁹⁻²¹

Figure 6. Evolution of cardiac troponin assays.¹⁹⁻²¹

Adapted from Apple FS, Collinson PO. *Clin Chem.* 2012;58:54-61; Thygesen K, et al. *Eur Heart J.* 2012;33:2252-2257
Al-Saleh A, et al. *CMAJ Open.* 2014;2:E199-207



The new generation of hs-cTn assays enables measurable values in a significant proportion of the reference population resulting in accurate assessment of the 99th percentile URL (CV ≤10%). The improved analytical sensitivity and precision of hs-cTn assays results in earlier detection of acute MI because of the reduction of the 'troponin-blind' interval with shortening of the window for serial testing to 3 hours or less.

3 Role of cardiac troponin in diagnosis of ACS

Interpretation and serial measurements

Troponin levels together with kinetic changes need to be integrated with clinical assessment and the ECG in the diagnosis of patients with suspected ACS (Figure 1). In its interpretation, cTn should be considered as a **quantitative marker of cardiac injury** with the likelihood of AMI increasing with the level of cTn (Figure 7).²²

Figure 7. cTn as a quantitative marker of cardiac injury.²²

Adapted from Wildi K, et al. *Clin Biochem*. 2015;48:218-222

Prediction of AMI		hsTn concentration	Likelihood of AMI (differential diagnosis)
Rule-in		Very high	Very large AMI myocarditis
	PPV >95%	...	Large AMI myocarditis, Takotsubo, PE, critical illness
	PPV 80%	High	Small AMI early large AMI, myocarditis, Takotsubo, PE, shock, heart failure, S aureus bacteremia (SAB), ...
	PPV 50% NPV 95%	...	Micro AMI early large AMI, myocarditis, Takotsubo, PE, shock, heart failure, hypertensive crisis, SAB, stable CAD, ...
----- 99 th percentile URL -----			
Rule-out	NPV 98%	...	Stable angina, heart failure, LVH, subclinical heart disease, ...
	NPV 99%	Very low	Healthy individuals

The lower the level of cTn, the lower the likelihood and the higher the negative predictive value (NPV) for the presence of AMI. In contrast, the higher the level, the higher the likelihood and the higher the positive predictive value (PPV) for the presence of AMI. Elevations above 5-fold the URL have a high (>90%) positive predictive value (PPV) for AMI. Elevations up to 3-fold the URL have only limited PPV (50-60%) for acute MI and may be associated with a broad spectrum of conditions.

Furthermore, it is important to consider the **difference between serial measurements (delta)** to differentiate chronic from acute conditions²²:

- Absolute changes in cTn provide a higher diagnostic accuracy for AMI as compared to relative changes. A caveat is that patients with AMI may show no or only a minimal change when troponin levels are measured around the peak of cTn release.
- The higher the absolute change, the higher the likelihood for AMI. However, there are several other acute cardiac conditions including tachyarrhythmias, myocarditis, hypertensive crisis, and Takotsubo cardiomyopathy that may also present with substantial cTn changes.

Algorithms for rule-in and rule-out of AMI

Based on a **decision limit** as well as criteria for a difference (**delta**) between serial measurement of cTn, **algorithms** have been designed and validated for rule-out and rule-in of AMI:

- Rule-out algorithms are optimized for **sensitivity** to minimize false-negatives.²³
- Rule-in algorithms are optimized for **specificity** to minimize false-positives.²⁴

Algorithms are based on the following **two approaches** (Table 5):

- 1. Fixed threshold based on the 99th percentile URL: decision within 3-6 hours from admission.**
This is the classical approach as recommended in the **Universal Definition of MI**.⁷ For the diagnosis of MI, serial measurement of cTn is needed for detection of a dynamic changing pattern of cTn with at least one value above the 99th percentile URL.
 - With **conventional cTn assays** a time window of **6 hours** or more from admission is needed.²⁵
 - Compared with conventional assays, **hs-cTn assays** allow earlier detection of acute MI, with shortening of the time window for serial testing to **3 hours**.^{5, 20}
- 2. Threshold NOT based on the 99th percentile URL: decision in 2 hours, or less, from admission.**

The introduction of hs-cTn assays, with accurate measurement below the 99th percentile URL and detection of small absolute changes within **1 or 2 hours**, has led to the development of rapid rule-out/rule-in strategies with substantial reduction in the time to decision.^{23,24}

Table 5. Algorithms for rule-out and rule-in of AMI for hs-cTn assays.^{23,24}

Adapted from Mueller C, et al. *Eur Heart J Acute Cardiovasc Care*. 2017;6:218-222; Möckel M, et al. *Eur Heart J Acute Cardiovasc Care*. 2017;6:212-217

	Decision limit	Serial measurement	Time for decision
Classical approach	99 th percentile URL	Yes	3 hours
Accelerated approach	Decision limits plus early absolute changes (cutoff levels are assay specific)	Yes (single sample rule-out on admission may also be used)	1 or 2 hours

Considerations for the use of cTn in diagnostic algorithms

Initial management of a patient with suspected ACS relies on clinical assessment **plus** 12-lead ECG for rapid detection of STEMI. If no immediate threat is detected, NSTEMI should be considered and the decision-making process is guided by triage algorithms that rely on measurement of cTn by **hs-cTn assays**.⁸

With hs-cTn assays the following triage algorithms are endorsed by ESC guidelines:⁸

- Standard algorithm (99th percentile URL): 0 h/3 h (Figure 8)
- Rapid algorithms: 0 h/1 h (Figure 9) or 0 h/2 h


Latest 2020 ESC Guidelines for the management of ACS in patients presenting without persistent ST-segment elevation recommend the use of an 0 h/2 h algorithm as an alternative to an 0 h/1 h algorithm.⁸

Because of lack of standardization, decision limits (thresholds) and delta values are assay specific. Consequently, algorithms are also assay-specific with a very low cut-off for rule-out, a high cut-off for rule-in and delta values based on the absolute change between two consecutive measurements to make further decisions on rule-out or rule-in. If hs-cTn is below the lower cutoff and remains low, the risk for 30-day adverse outcome is low and AMI is ruled out. If the cTn concentration is above the high threshold, the probability for AMI is high and the patient should undergo CCTA within 48 hours (see flow chart in Figure 10).

Algorithms should always be used in conjunction with all available clinical information.⁸ Special attention should be paid to time of chest pain onset as it impacts interpretation of the results and the need for further serial measurements. For very early presenters (<1 h from chest pain onset) an additional sample will be needed at 3 hours.⁸ Late increases of cTn may occur in 1% of patients. Therefore, repeat testing at later time points (3-6 hours) is advised if first measurements are not conclusive and clinical suspicion remains high.⁸

Figure 8. ESC 0 h/3 h algorithm.⁵

Adapted from Roffi M, et al. *Eur Heart J*. 2016;37:267-315




GRACE = Global Registry of Acute Coronary Events score; hs-cTn = sensitivity cardiac troponin; URL = upper reference limit, 99th percentile of healthy controls.

* Δ change, dependent on assay. Highly abnormal hsTn defines values beyond 5-fold the upper reference limit.

Valid for hs-cTn assays in patients presenting with suspected NSTEMI to the ED. Cut-off is based on the assay-specific 99th percentile URL. An increase of >50% of the URL is considered a relevant delta change.²⁰

Figure 9. Rapid 0 h/1 h algorithm.⁸


Adapted from Collet JP, et al. *Eur Heart J*. 2021;42:1289-1367



Valid for hs-cTn assays in patients presenting with suspected NSTEMI to the ED. Cut-off levels and delta values are assay-specific. Alternative rapid 0 h/2 h algorithms also exist.

Figure 10. Proposed flow chart for rapid algorithms.

Adapted from Collinson P. *Clin Biochem.* 2021;S0009-9120(21)00047-3



Follow-up medical decisions for patients categorized as high risk (rule-in), intermediate risk and low risk (rule-out).

ACS, acute coronary syndrome; cTn, cardiac troponin; ECG, electrocardiogram; LoD, limit of detection; NSTEMI, non-ST-segment elevation myocardial infarction; STEMI, ST-segment-elevation myocardial infarction; URL, upper reference limit

4 Guide to interpretation

Pre-analytical and analytical considerations

➔ Sample type

Serum, plasma, or anticoagulated **whole blood** (point-of-care tests) are acceptable specimens for the analysis of cTn, but different specimen types should not be used interchangeably.²⁶ Depending on the assay, differences may exist between serum and plasma samples.²⁵ Also, there may be an influence of the type of anticoagulant that is used for sample collection with lower values in EDTA- vs. lithium heparin plasma.²⁵

Rapid processing of samples is important and **results** should be obtained **within 60 minutes of blood collection**.⁸ Therefore, whole blood and plasma are the specimens of choice for the emergency laboratory to avoid delays due to serum preparation.²⁵

➔ Precision

More sensitive assays provide better precision at low values and facilitate the detection of changing values within 1 or 2 hours of admission. Key criteria for a hs-cTn assay are²⁷:

- Analytical imprecision (% CV) at its established 99th percentile URL $\leq 10\%$.
- Ability to quantify troponin levels in the majority (>50%) of healthy subjects.

➔ Standardization

In contrast to cTnT, with only one supplier, there are many cTnI assays commercially available. These assays use different calibrators and have a variable immunoreactivity to the various circulating forms of cTnI. Consequently, the standardization of cTnI assays remains a huge challenge.²⁸ Therefore, despite correlation between methods, reported values of cTnI results differ across measurement platforms and each assay will have its own assay specific clinical decision limits.

Expected values

➔ Impact of age, gender and ethnicity

With hs-cTn assays a difference between men and women becomes apparent with higher values in men compared with women.²⁹ This is biologically plausible because of the larger heart mass in men. Furthermore, in adults cTn increases with age in both men and women.²⁹ The reverse trend is seen in healthy children, with highest values in the first month of life and a gradual decline until adolescence.³⁰ Concerning ethnicity, higher values have been reported in African-Americans compared with Caucasians.²⁹

Laboratory practice guidelines recommend that manufacturers report sex-specific 99th percentile URL's for their hs-cTn assay.²⁶ Specific URLs by age or by ethnicity are currently not recommended.²⁶ Despite clear sex differences in URL's, the clinical superiority of sex-specific URL's have not been definitely proven.³¹

➔ Causes of elevated cTn

It is important to consider cTn as a **quantitative marker of cardiac injury**: the higher the level, the higher the likelihood of AMI (**Figure 7**).

Elevated cTn (i.e. above the 99th percentile URL) is not specific for AMI and can be observed in numerous conditions.

First, cTn is not only elevated in the classical **type 1 MI** (atherothrombosis) but also in **type 2 MI** due to a myocardial oxygen supply-demand mismatch associated with a variety of underlying disorders (**Table 2**).

Besides type 1 and type 2 MI, there are a number of other cardiac and non-cardiac conditions where elevated cTn levels can occur (**Table 2**).

Elevated cTn in critical care is indicative of myocardial injury and an adverse prognosis.³²

This has also been reported for the most recent respiratory pathogen SARS-CoV-2 (COVID-19).³³

FREQUENTLY ASKED QUESTIONS

1 What is the difference between cTnI and cTnT?

Guidelines do not differentiate between cTnI and cTnT as a tool in the diagnosis and prognosis of suspected ACS.^{5,7,8} In a direct comparison both hs-cTnI and hs-cTnT provided overall similar high diagnostic and prognostic accuracy, albeit that cTnI showed a somewhat higher diagnostic accuracy for early presenters and cTnT for late presenters, and also cTnT showed a higher prognostic accuracy for all-cause mortality.³⁴

The following are **notable differences between cTnI and cTnT**:

- After AMI, cTnT increases tend to continue for longer than cTnI increases.³⁵
- Renal dysfunction leads more frequently to elevations in cTnT than cTnI.³⁶
- Elevated cTnT, and not cTnI, may occur in chronic skeletal muscle damage (e.g. muscular dystrophies or chronic myositis).³⁶
- cTnT, but not cTnI, shows a diurnal rhythm characterized by gradually decreasing concentrations throughout daytime, rising concentrations during night-time, with peak concentrations in the morning.³⁷
- There is only one supplier of cTnT assays but multiple suppliers of cTnI assays, each with their own monoclonal antibodies and calibrators. This variability in source materials, combined with a lack of standardization, results in a large heterogeneity and poor concordance between cTnI assays.³⁸ Consequently, for cTnI, results are not transferable from one assay to another.

2 Can false-negative or false-positive cTn results be expected?

A false-positive result is a test result above the decision threshold in the absence of disease, whereas a false-negative result is a test result below the decision threshold in the presence of disease. However, it is important to realize that although cTn is highly specific for cardiac damage it has only moderate specificity for AMI and can be elevated in many other clinical settings (**Table 2**). Here, the term false-positive is misleading because in these clinical settings an elevated cTn test result, as a marker of cardiac damage, is associated with adverse outcomes.³⁹

The importance of a significant change in cTn levels upon serial testing is a key consideration in distinguishing AMI from other causes of cardiac damage. Furthermore, serial testing is also important to avoid false-

negatives which may occur in very early presenters or due to late increases in cTn. Therefore, repeat testing at later time points is always needed if first measurements are not conclusive and clinical suspicion remains high.⁸

Analytical interferences are a rare cause of false results but should be considered when test results and clinical presentation are discordant.³⁶ Most often these issues cause increased cTn values that do not change acutely over time. Possible reasons for erroneous results include the following:³⁶

- Analytical interferences due to pronounced hemolysis, hyperlipidemia or hyperbilirubinemia.
- Non-reproducible increases in concentrations, referred to as 'flyers' or 'outliers'.
- Presence of heterophilic or human anti-mouse antibodies.
- Presence of auto-antibodies against cTnI or cTnT. These antibodies can directly interfere with the assay or their binding to cTn can result in high molecular weight complexes (macrotroponins) that are cleared more slowly from the circulation, resulting in persistently elevated cTn concentrations.

In the case of erroneous and discordant results, it is important to consult the lab to investigate and eliminate potential analytical interferences.

3 How is the 99th percentile URL for cTn assays determined?

The 99th percentile URL need to be determined for both men and women which requires at least 300 normal adult individuals of each sex with a diverse distribution of ages from 18 to over 70 years.^{27,40} Apart from age and gender, cTn levels may be influenced by underlying factors with impact on cardiac structure and function.⁴¹ A critical element, therefore, is the definition of 'normal', because studies have shown that a more rigorous selection based on eliminating potential comorbidities will lower the URL.^{40,42}

Currently, there is no uniform consensus on how to select 'cardiac healthy' individuals, but experts have proposed the following **exclusion criteria**:^{26,40}

- History of cardiovascular disease and medication usage.
- Diabetes (HbA1c $\geq 6.5\%$).
- Renal dysfunction (eGFR < 60 mL / min / 1.73 m²).
- Elevated NT-proBNP as a surrogate marker of cardiac dysfunction (>125 pg/mL for age < 75 years and >450 pg/mL for age ≥ 75 years).

4 What are the clinical implications of high-sensitivity cTn assays?

The clinical implications of hs-cTn assays are summarized in **Table 6**. Compared with conventional cTn assays, the new generation of hs-cTn assays have **improved diagnostic and prognostic properties** with the following **impact on clinical management**:

→ **Faster time to decision**

The improved diagnosis translates into a reduction of the 'troponin-blind' interval leading to earlier detection of acute MI. This particularly enables earlier rule-out because of the higher NPV for acute MI at presentation.²¹ The introduction of hs-cTn was reported to be associated with an improved rule-out process with less need for stress testing, faster time to discharge and a reduction in costs.⁴³

→ **Change in case mix**

Introduction of hs-cTn assays results in about a 4% increase in the incidence of AMI with a corresponding decrease in the diagnosis of UA.^{8, 44, 45} Furthermore, there is a 2-fold increase in the detection of type 2 MI.⁸

→ **Better risk stratification**

Introduction of hs-cTn assays has been shown to identify more patients at risk for mortality who may be suitable for beneficial cardiovascular therapies.^{44, 46, 47} In Sweden, this resulted in an increased use of coronary angiographies and revascularizations with no impact on survival but a reduced risk of reinfarction.⁴⁵

→ **Better patient outcomes**

In patients with suspected ACS, implementation of a sensitive cTn assay increased the diagnosis of MI and identified patients at high risk of recurrent MI and death. Lowering the diagnostic threshold of cTn was associated with major reductions in morbidity and mortality.⁴⁸

Table 6. Clinical implications of high-sensitivity cardiac troponin assays.⁸

Adapted from Collet JP, et al. *Eur Heart J.* 2021;42:1289-1367

Compared with standard cTn assays, hs-cTn assays

- Have higher NPV for AMI.
- Reduce the 'troponin blind' interval leading to earlier detection of AMI.
- Result in ~4% absolute and ~20% relative increases in detection of type 1 MI and a corresponding decrease in the diagnosis of UA.
- Are associated with a 2-fold increase in the detection of type 2 MI.

Levels of hs-cTn should be interpreted as quantitative markers of cardiomyocyte damage (i.e. the higher the level, the greater the likelihood of MI):

- Elevations beyond 5-fold the URL have high (>90%) PPV for acute type 1 MI.
- Elevations beyond 3-fold the URL have only limited (50-60%) PPV for AMI and may be associated with a broad spectrum of conditions.
- It is common to detect circulating levels of cTn in healthy individuals

Rising and/or falling cTn levels differentiate acute (as in MI) from chronic cardiomyocyte damage (the more pronounced the change, the higher the likelihood of AMI)

5 Is there a role for cTn in the distinction between type 1 and type 2 MI?

Following the introduction of hs-cTn assays the incidence of type 2 MI and ischemic myocardial injury is projected to increase.⁴⁹ The distinction between type 1 and type 2 MI and ischemic myocardial injury remains a challenge but is critically important because of differences in therapy.⁴⁹ The differential diagnosis between type 1 and 2 MI can only be done with coronary angiography, but such invasive approach is not warranted or can even be harmful in type 2 MI patients. Therefore, there is a clinical need for a rapid non-invasive differentiation of both MI types.

In general, cTn values are higher in type 1 than in type 2 MI.⁵⁰ Yet, no cut-off has been defined to accurately discriminate between both types and also the use of delta values is not helpful.⁵¹ Compared with type 1 MI, patients with type 2 MI are more often female and are more often not presenting with typical radiating chest pain.⁵¹ These observations have led to the development of a simple binary score to predict the probability of type 2 MI, with highest probability (72%) when all 3 predictors are positive: female sex, absence of radiating chest pain and hs-cTnI ≤40.8 ng/L.⁵¹ This score still requires external validation, including appropriateness of this cut-off with other hs-cTnI assays.

6 How do I select the most optimal rule-out/ rule-in algorithm?

Several elements are important to consider in this selection process. In the ED setting of patients with chest pain and suspected ACS the most valuable application is rapid and safe rule-out. This requires an informed balance between safety (NPV for the presence of AMI) and efficacy (proportion triaged towards rule-out). Furthermore, it is important to consider the user friendliness of the choice of algorithm and the impact on patient flow and staff workload, hence it is dependent on the process of the particular hospital.

The latest 2020 ESC Guidelines for the management of ACS in patients presenting without persistent ST-segment elevation clearly recommend the use of a 0 h/2 h algorithm as an alternative to a validated 0 h/1 h algorithm.⁸ The accelerated algorithm implemented on a particular site will depend on the test available at the central laboratory (hs-cTnI or hs-cTnT). It will also be influenced by the ability of the laboratory to return results within 45 minutes of blood draw.

Both accelerated (0 h/1 h or 0 h/2 h) and 0 h/3 h algorithms are recommended by the ESC for early rule-out of AMI.⁸ These are completely different protocols. The accelerated algorithms are based on hs-cTn levels at

presentation and absolute changes within the first 1 or 2 hours. They are based on the ability to predict that the patient will rule-in or rule-out for MI according to the Universal Definition of MI on further testing. This is discussed in more detail below. The 0 h/3 h algorithm uses the fixed 99th percentile URL threshold in conjunction with clinical criteria (GRACE score) and is based directly on the Universal Definition of MI. From a direct comparison of both approaches it was concluded that the accelerated algorithms are superior because they more favorably combine safety with efficacy.⁵² Advantages of the accelerated algorithms are that they allow earlier clinical decision making and do not require a specific risk score which increases its feasibility.

The 0 h/1 h approach has been questioned on the basis of patient flow and ED and laboratory logistics related to sampling at this short time interval. Consequently, a 2-hour protocol may be more suitable for some institutions.⁵³

7 What are the limitations of rapid rule-out/ rule-in algorithms?

As a critical principle it is important to emphasize that AMI remains a clinical diagnosis and that the results of cTn testing need to be interpreted in conjunction with all other available information such as 12-lead ECG, patient history and examination and other diagnostic investigations. Rapid cTn testing algorithms are very useful to triage about 70% of suspected patients towards either AMI rule-out or AMI rule-in. This is not synonymous with a final diagnosis for presence or absence of AMI, but only informs about the probability of having AMI which aids in making decisions on the next course of action in the clinical management of the patient.

➔ AMI rule-out

In this case, the decision to discharge home from the ED must not be solely based on cTn testing. Repeat cTn testing at later time points may be required if first measurements are not conclusive and clinical suspicion for ACS remains high, for example in the case of early presenters (<2 hours from onset of chest pain).^{5,8} AMI rule-out does not mean ACS rule-out and, therefore, further investigations such as functional cardiac stress testing may be needed to rule-out UA. Also, in the case of chest pain, it is important to be aware of other life-threatening disorders such as AAD or PE.

In the context of AMI rule-out, a **clinical risk score** may provide incremental value for the selection of patients who are also candidates for **rapid discharge** from the ED and outpatient management.⁵⁴ For this reason, the **GRACE score** is already integrated in the 3-hour ESC algorithm (**Figure 8**).⁵ Compared with a 2-hour rule-out protocol, a 2-hour protocol that combined cTn testing with the **TIMI score** was superior for selection of patients for early discharge from the ED.⁵⁴

The **HEART score** outperformed both TIMI and GRACE in its capacity to predict major adverse cardiac events (MACE) after 30 days from presentation to the ED.⁵⁵ The HEART score integrates cTn testing with ECG, age and clinical predictors.⁵⁶ A single cTn at presentation and a HEART score of ≤ 3 demonstrated NPV of $\geq 99.5\%$ for 30-day MACE, which suggests that this a threshold for safe discharge.⁵⁵

➔ AMI rule-in

In this case the PPV is only about 75-80%. This is enough as a working diagnosis and is helpful for the physician to make a decision for coronary angiography and start of certain medical therapies. Conditions like Takotsubo cardiomyopathy and myocarditis can also be associated with substantial elevations and changes in cTn (**Figure 7**). To differentiate AMI from these conditions early coronary angiography is the general clinical consequence of triage towards AMI rule-in.

➔ Observational zone

Diagnostic uncertainty remains in approximately 25-30% of suspected patients.⁵⁷ Patients in the observational zone are typically elderly men with pre-existing CAD and have a similar high long-term mortality as patients triaged towards AMI rule-in.⁵⁸ Additional diagnostic investigations are needed in this group, including cTn testing at later time points (3-6 hours), functional stress imaging and coronary angiography in order to discriminate ACS (21% UA, 15% AMI) from non-cardiac disease (38%) and non-coronary cardiac disease (24%).⁵⁸

8 Should early presenters be managed differently?

Early presenters are patients with suspected ACS that are seen in the ED <2 hours from chest pain onset. High-sensitivity assays detect changes in troponin earlier and therefore improve the diagnostic accuracy for MI in these patients.

In early presenters, better diagnostic accuracy has been reported for cTnI as compared to cTnT.

The 0/3h algorithm allows longer time for clinical assessment and later sequential measurements, and will safely rule out ACS patients among early presenters.

The ESC 0 h/1h and 0 h/2h algorithms can also be used in these patients as safety (negative predictive value and sensitivity) is maintained for these strategies.

However, due to the kinetics of cTn release following MI, it is recommended that sequential measurements of hs-cTn are used in patients that present very early (<1h) after chest pain onset, especially if the first result is inconclusive and there is clinical suspicion of ACS. Although changes within the reference interval can be detected within 1-2 hours of MI, values above the 99th percentile can be reliably detected only 3 hours or more from symptom onset.

9 Which patients pose a challenge for the diagnosis of AMI by cTn?

Special attention is required for suspected patients presenting with renal dysfunction (eGFR <60 mL / min / 1.73 m²) or symptoms other than typical chest pain such as dyspnea.

➔ **Renal dysfunction** may occur in 16% of patients presenting with symptoms of suspected AMI.⁵⁹ These patients more frequently present with atypical symptoms and the proportion of diagnosed AMI is 2-fold higher than in suspected patients without renal function.⁵⁹ These patients pose a diagnostic challenge because, even in the absence of AMI, ECG abnormalities (ST-segment depression/T-wave inversion) and elevations of cTn above the 99th percentile URL are more frequently observed.⁵⁹ Yet, the diagnostic utility of hs-cTn assays is retained in patients with renal dysfunction as long as test specific optimized cut-off levels are used (2 to 3-fold higher than 99th percentile URL).⁵⁹

➔ **Dyspnea** may occur in 17-49% of patients with suspected ACS but the prevalence of diagnosed ACS is much lower than in patients presenting with typical chest pain.⁶⁰ This high-risk cohort poses a challenge for the diagnosis of AMI because of the low pre-test probability and the lower specificity of cTn.⁶⁰ Higher cut-offs at presentation and serial measurements over 6 hours may help to identify NSTEMI in this population.⁶⁰

10 What is the prognostic utility of cTn?

It is recommended to use established risk scores (e.g. GRACE) for prognosis estimation.⁸ These scores integrate hs-cTn test results with clinical factors.⁸ Invasive coronary angiography has a central role in the management of patients with NSTEMI-ACS and the timing of this intervention is guided by risk stratification. An early invasive strategy, defined as coronary angiography within 24 hours of hospital admission, is required in the presence of high-risk criteria, including a dynamic rise or fall in cTn compatible with MI.⁸

Irrespective of the presence of AMI, any elevated cTn is predictive of an adverse outcome.³⁹ The HEART score, which integrates cTn and hs-cTn testing with ECG, age and clinical predictors, is a useful tool to select MI rule-out patients for safe discharge.⁵⁵

11 Is there a role for cTn testing in other diseases?

Information on an elevated cTn (>99th percentile URL) may be helpful in the management of the following cardiovascular diseases:

➔ Heart failure⁶¹

Elevated cTn is frequently observed in both acute and chronic HF due to cardiac or non-cardiac mechanisms. In the case of acute HF, cTn should be measured to confirm or exclude MI as the precipitating cause. In this setting, an elevated cTn is associated with a more decompensated profile and a higher risk for an adverse outcome. The 10% of acute HF patients with cTn below the 99th percentile URL are at very low risk for a short-term adverse event and may be eligible for early discharge.⁶² An elevated cTn has independent prognostic value in chronic HF, but data are lacking about specific therapeutic interventions in such cases.

➔ Pulmonary embolism⁶³

Approximately 50% of acute PE patients show elevated cTn values on admission and these patients have a worse prognosis. Decisions on clinical management of PE in the acute phase are guided by the assessment of the patient's early death risk. This risk assessment integrates clinical parameters, assessment of right ventricular dysfunction (RVD) by echocardiography and measurement of cardiac biomarkers (cTn or natriuretic peptides). Absence of an elevated cTn has a high NPV (>98%) for a complicated course. Patients stratified as low risk based on a low clinical score, absence of RVD and normal cTn may be potential candidates for early discharge and home treatment.

➔ Atrial fibrillation^{64, 65}

Currently, a clinical score is recommended for stroke risk prediction in patients with AF in order to select patients for oral anticoagulant therapy. In AF patients on oral anticoagulation it is important to identify modifiable risk factors for major bleeding. Higher concentrations of hs-cTn in AF patients are independently associated with increased risk for stroke, cardiac mortality or bleeding. Because of this additional prognostic information in AF patients, a biomarker-based risk score for stroke as well as for bleeding may prove helpful for better risk stratification.

REFERENCES

- Roth GA, Mensah GA, Johnson CO, et al. GBD-NHLBI-JACC Global Burden of Cardiovascular Diseases Writing Group. **Global burden of cardiovascular diseases and risk factors, 1990-2019: update from the GBD 2019 Study.** *J Am Coll Cardiol.* 2020;76:2982-3021.
- Timmis A, Townsend N, Gale CP, et al. **European Society of Cardiology: Cardiovascular Disease Statistics 2019.** *Eur Heart J.* 2020;41:12-85.
- GBD 2019 Diseases and Injuries Collaborators. **Global burden of 369 diseases and injuries in 204 countries and territories, 1990-2019: a systematic analysis for the Global Burden of Disease Study 2019.** *Lancet.* 2020;396:1204-1222.
- Rui P, Kang K. **National Hospital Ambulatory Medical Care Survey: 2017 emergency department summary tables.** National Center for Health Statistics. Available from: https://www.cdc.gov/nchs/data/nhamcs/web_tables/2017_ed_web_tables-508.pdf
- Roffi M, Patrono C, Collet JP, et al.; ESC Scientific Document Group. **2015 ESC Guidelines for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation: Task Force for the Management of Acute Coronary Syndromes in Patients Presenting without Persistent ST-Segment Elevation of the European Society of Cardiology (ESC).** *Eur Heart J.* 2016;37:267-315.
- Ibanez B, James S, Agewall S, et al.; ESC Scientific Document Group. **2017 ESC Guidelines for the management of acute myocardial infarction in patients presenting with ST-segment elevation: The Task Force for the management of acute myocardial infarction in patients presenting with ST-segment elevation of the European Society of Cardiology (ESC).** *Eur Heart J.* 2018;39:119-177.
- Thygesen K, Alpert JS, Jaffe AS, et al.; Executive Group on behalf of the Joint European Society of Cardiology (ESC)/American College of Cardiology (ACC)/American Heart Association (AHA)/World Heart Federation (WHF) Task Force for the Universal Definition of Myocardial Infarction. **Fourth Universal Definition of Myocardial Infarction (2018).** *Circulation.* 2018;138:e618-e651.
- Collet JP, Thiele H, Barbato E, et al.; ESC Scientific Document Group. **2020 ESC Guidelines for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation.** *Eur Heart J.* 2021;42:1289-1367.
- Sandoval Y, Thygesen K. **Myocardial infarction type 2 and myocardial injury.** *Clin Chem.* 2017;63:101-107.
- Chapman AR, Shah ASV, Lee KK, et al. **Long-term outcomes in patients with type 2 myocardial infarction and myocardial injury.** *Circulation.* 2018;137:1236-1245.
- Partida RA, Libby P, Crea F, Jang IK. **Plaque erosion: a new in vivo diagnosis and a potential major shift in the management of patients with acute coronary syndromes.** *Eur Heart J.* 2018;39:2070-2076.
- Shah ASV, Sandoval Y, Noaman A, et al. **Patient selection for high sensitivity cardiac troponin testing and diagnosis of myocardial infarction: prospective cohort study.** *BMJ.* 2017; 359:j4788.
- Swap CJ, Nagurney JT. **Value and limitations of chest pain history in the evaluation of patients with suspected acute coronary syndromes.** *JAMA.* 2005;294:2623-2629.
- Fanaroff AC, Rymer JA, Goldstein SA, et al. **Does this patient with chest pain have acute coronary syndrome?: The rational clinical examination systematic review.** *JAMA.* 2015;314:1955-1965.
- Parmacek MS, Solaro RJ. **Biology of the troponin complex in cardiac myocytes.** *Prog Cardiovasc Dis.* 2004;47:159-176.
- Mair J, Lindahl B, Hammarsten O, et al.; European Society of Cardiology (ESC) Study Group on Biomarkers in Cardiology of the Acute Cardiovascular Care Association (ACCA). **How is cardiac troponin released from injured myocardium?** *Eur Heart J Acute Cardiovasc Care.* 2018;7:553-560.
- Bates KJ, Hall EM, Fahie-Wilson MN, et al. **Circulating immunoreactive cardiac troponin forms determined by gel filtration chromatography after acute myocardial infarction.** *Clin Chem.* 2010;56:952-958.
- Apple FS, Sandoval Y, Jaffe AS, Ordonez-Llanos J; IFCC Task Force on Clinical Applications of Cardiac Bio-Markers. **Cardiac troponin assays: guide to understanding analytical characteristics and their impact on clinical care.** *Clin Chem.* 2017;63:73-81.
- Apple FS, Collinson PO; IFCC Task Force on Clinical Applications of Cardiac Biomarkers. **Analytical characteristics of high-sensitivity cardiac troponin assays.** *Clin Chem.* 2012;58:54-61.
- Thygesen K, Mair J, Giannitsis E, et al.; Study Group on Biomarkers in Cardiology of ESC Working Group on Acute Cardiac Care. **How to use high-sensitivity cardiac troponins in acute cardiac care.** *Eur Heart J.* 2012;33:2252-2257.
- Al-Saleh A, Alazzonei A, Al Shalash S, et al. **Performance of the high-sensitivity troponin assay in diagnosing acute myocardial infarction: systematic review and meta-analysis.** *CMAJ Open.* 2014;2:E199-207.
- Wildi K, Twerenbold R, Mueller C. **How acute changes in cardiac troponin concentrations help to handle the challenges posed by troponin elevations in non-ACS-patients.** *Clin Biochem.* 2015;48:218-222.
- Mueller C, Giannitsis E, Möckel M, et al.; Biomarker Study Group of the ESC Acute Cardiovascular Care Association. **Rapid rule out of acute myocardial infarction: novel biomarker-based strategies.** *Eur Heart J Acute Cardiovasc Care.* 2017;6:218-222.
- Möckel M, Giannitsis E, Mueller C, et al.; Biomarker Study Group of the European Society of Cardiology Acute Cardiovascular Care Association. **Editor's Choice - Rule-in of acute myocardial infarction: Focus on troponin.** *Eur Heart J Acute Cardiovasc Care.* 2017;6:212-217.
- Thygesen K, Mair J, Katus H, et al.; Study Group on Biomarkers in Cardiology of the ESC Working Group on Acute Cardiac Care. **Recommendations for the use of cardiac troponin measurement in acute cardiac care.** *Eur Heart J.* 2010;31:2197-2204.
- Wu AHB, Christenson RH, Greene DN, et al. **Clinical Laboratory practice recommendations for the use of cardiac troponin in acute coronary syndrome: Expert opinion from the Academy of the American Association for Clinical Chemistry and the Task Force on Clinical Applications of Cardiac Bio-Markers of the International Federation of Clinical Chemistry and Laboratory Medicine.** *Clin Chem.* 2018;64:645-655.
- Apple FS, Jaffe AS, Collinson P, et al.; International Federation of Clinical Chemistry (IFCC) Task Force on Clinical Applications of Cardiac Bio-Markers. **IFCC educational materials on selected analytical and clinical applications of high sensitivity cardiac troponin assays.** *Clin Biochem.* 2015;48:201-203.
- Christenson RH, Bunk DM, Schimmel H, Tate JR; IFCC Working Group on Standardization of Troponin I. **Point: Put simply, standardization of cardiac troponin I is complicated.** *Clin Chem.* 2012;58:165-168.
- Gore MO, Seliger SL, Defilippi CR, et al. **Age- and sex-dependent upper reference limits for the high-sensitivity cardiac troponin T assay.** *J Am Coll Cardiol.* 2014;63:1441-1448.

30. Caselli C, Cangemi G, Masotti S, et al. **Plasma cardiac troponin I concentrations in healthy neonates, children and adolescents measured with a high sensitive immunoassay method: High sensitive troponin I in pediatric age.** *Clin Chim Acta.* 2016;458:68-71.
31. Eggers KM, Lindahl B. **Impact of sex on cardiac troponin concentrations - A critical appraisal.** *Clin Chem.* 2017;63:1457-1464.
32. Reynolds T, Ceconi M, Collinson P, et al. **Raised serum cardiac troponin I concentrations predict hospital mortality in intensive care unit patients.** *Br J Anaesth.* 2012;109:219-224.
33. Li X, Pan X, Li Y, et al. **Cardiac injury associated with severe disease or ICU admission and death in hospitalized patients with COVID-19: a meta-analysis and systematic review.** *Crit Care.* 2020;24:468.
34. Rubini Gimenez M, Twerenbold R, Reichlin T, et al. **Direct comparison of high-sensitivity-cardiac troponin I vs. T for the early diagnosis of acute myocardial infarction.** *Eur Heart J.* 2014;35:2303-2311.
35. Mair J, Thome-Kromer B, Wagner I, et al. **Concentration time courses of troponin and myosin subunits after acute myocardial infarction.** *Coron Artery Dis.* 1994; 5: 865-872.
36. Mair J, Lindahl B, Müller C, et al. **What to do when you question cardiac troponin values.** *Eur Heart J Acute Cardiovasc Care.* 2018;7:577-586.
37. Klinkenberg LJ, Wildi K, van der Linden N, et al. **Diurnal rhythm of cardiac troponin: consequences for the diagnosis of acute myocardial infarction.** *Clin Chem.* 2016;62:1602-1611.
38. Chenevier-Gobeaux C, Deweerdt L, Cantero AV, et al.; Troponins Working Group of the Société Française de Biologie Clinique (SFBC), in collaboration with the Société Française de Médecine d'Urgence (SFMU) and the Société Française de Cardiologie (SFC). **Multi-centre evaluation of recent troponin assays for the diagnosis of NSTEMI.** *Pract Lab Med.* 2018;11:23-32.
39. Giannitsis E, Katus HA. **Cardiac troponin level elevations not related to acute coronary syndromes.** *Nat Rev Cardiol.* 2013;10:623-634.
40. Sandoval Y, Apple FS. **The global need to define normality: the 99th percentile value of cardiac troponin.** *Clin Chem.* 2014;60:455-462.
41. de Lemos JA, Grundy SM. **Low levels of circulating troponin as an intermediate phenotype in the pathway to heart failure.** *J Am Coll Cardiol.* 2012;59:490-492.
42. Collinson PO, Heung YM, Gaze D, et al. **Influence of population selection on the 99th percentile reference value for cardiac troponin assays.** *Clin Chem.* 2012;58:219-225.
43. Twerenbold R, Jaeger C, Rubini Gimenez M, et al. **Impact of high-sensitivity cardiac troponin on use of coronary angiography, cardiac stress testing, and time to discharge in suspected acute myocardial infarction.** *Eur Heart J.* 2016;37:3324-3332.
44. Reichlin T, Twerenbold R, Reiter M, et al. **Introduction of high-sensitivity troponin assays: impact on myocardial infarction incidence and prognosis.** *Am J Med.* 2012;125:1205-1213.e1.
45. Odqvist M, Andersson PO, Tygesen H, et al. **High-sensitivity troponins and outcomes after myocardial infarction.** *J Am Coll Cardiol.* 2018;71:2616-2624.
46. Melki D, Lugnegård J, Alfreðsson J, et al. **Implications of introducing high-sensitivity cardiac troponin T into clinical practice: data from the SWEDEHEART Registry.** *J Am Coll Cardiol.* 2015;65:1655-1664.
47. Eggers KM, Lindahl B, Melki D, Jernberg T. **Consequences of implementing a cardiac troponin assay with improved sensitivity at Swedish coronary care units: an analysis from the SWEDEHEART registry.** *Eur Heart J.* 2016;37:2417-24.
48. Mills NL, Churchhouse AMD, Lee KK, et al. **Implementation of a sensitive troponin I assay and risk of recurrent myocardial infarction and death in patients with suspected acute coronary syndrome.** *JAMA.* 2011;305:1210-1216.
49. McCarthy CP, Vaduganathan M, Januzzi JL Jr. **Type 2 myocardial infarction-diagnosis, prognosis, and treatment.** *JAMA.* 2018;320:433-434.
50. Lippi G, Sanchis-Gomar F, Cervellin G. **Cardiac troponins and mortality in type 1 and 2 myocardial infarction.** *Clin Chem Lab Med.* 2017;55:181-188.
51. Sandoval Y, Thorsden SE, Smith SW, et al. **Cardiac troponin changes to distinguish type 1 and type 2 myocardial infarction and 180-day mortality risk.** *Eur Heart J Acute Cardiovasc Care.* 2014;3:317-325.
52. Badertscher P, Boeddinghaus J, Twerenbold R, et al.; APACE Investigators. **Direct comparison of the 0/1h and 0/3h algorithms for early rule-out of acute myocardial infarction.** *Circulation.* 2018;137:2536-2538.
53. Boeddinghaus J, Reichlin T, Cullen L, et al. **Two-hour algorithm for triage toward rule-out and rule-in of acute myocardial infarction by use of high-sensitivity cardiac troponin I.** *Clin Chem.* 2016;62:494-504.
54. Wildi K, Cullen L, Twerenbold R, et al. **Direct comparison of 2 rule-out strategies for acute myocardial infarction: 2-h accelerated diagnostic protocol vs 2-h algorithm.** *Clin Chem.* 2017;63:1227-1236.
55. Reaney PDW, Elliot HI, Noman A, Cooper JG. **Risk stratifying chest pain patients in the emergency department using HEART, GRACE and TIMI scores, with a single contemporary troponin result, to predict major adverse cardiac events.** *Emerg Med J.* 2018;35:420-427.
56. www.heartscore.nl. Accessed on Feb 03 2021.
57. Lindahl B, Jernberg T, Badertscher P, et al. **An algorithm for rule-in and rule-out of acute myocardial infarction using a novel troponin I assay.** *Heart.* 2017;103:125-131.
58. Nestelberger T, Wildi K, Boeddinghaus J, et al. **Characterization of the observe zone of the ESC 2015 high-sensitivity cardiac troponin 0h/1h-algorithm for the early diagnosis of acute myocardial infarction.** *Int J Cardiol.* 2016;207:238-245.
59. Twerenbold R, Wildi K, Jaeger C, et al. **Optimal cutoff levels of more sensitive cardiac troponin assays for the early diagnosis of myocardial infarction in patients with renal dysfunction.** *Circulation.* 2015;131:2041-2050.
60. Biener M, Mueller M, Vafaie M, et al. **Impact of leading presenting symptoms on the diagnostic performance of high-sensitivity cardiac troponin T and on outcomes in patients with suspected acute coronary syndrome.** *Clin Chem.* 2015;61:744-751.
61. Januzzi JL Jr, Filippatos G, Nieminen M, Gheorghiade M. **Troponin elevation in patients with heart failure: on behalf of the third Universal Definition of Myocardial Infarction Global Task Force: Heart Failure Section.** *Eur Heart J.* 2012;33:2265-2271.
62. Jaffe AS, Morrow DA, Scirica BM. **High-sensitivity troponin in the triage of acute decompensated heart failure.** *JACC Heart Fail.* 2016;4:600-602.
63. Konstantinides SV, Torbicki A, Agnelli G, et al.; Task Force for the Diagnosis and Management of Acute Pulmonary Embolism of the European Society of Cardiology (ESC). **2014 ESC guidelines on the diagnosis and management of acute pulmonary embolism.** *Eur Heart J.* 2014;35:3033-3069, 3069a-3069k.
64. Kirchhof P, Benussi S, Kotecha D, et al.; ESC Scientific Document Group. **2016 ESC Guidelines for the management of atrial fibrillation developed in collaboration with EACTS.** *Eur Heart J.* 2016;37:2893-2962.
65. Hijazi Z, Oldgren J, Siegbahn A, Wallentin L. **Application of biomarkers for risk stratification in patients with atrial fibrillation.** *Clin Chem.* 2017;63:152-164.

LIST OF ABBREVIATIONS & ACRONYMS

AAD	acute aortic dissection
ACS	acute coronary syndromes
AF	atrial fibrillation
AMI	acute myocardial infarction
CABG	coronary artery bypass grafting
CAD	coronary artery disease
CCTA	coronary computed tomography angiography
CCU	coronary care unit
cTn	cardiac troponin
CV	coefficient of variation
ECG	electrocardiogram
ED	emergency department
EDTA	ethylene diamine tetra-acetic acid
eGFR	estimated glomerular filtration rate
ESC	European Society of Cardiology
GRACE score	Global Registry of Acute Coronary Events; the GRACE risk score has been developed to assess the risk of mortality in patients with ACS www.mdcalc.com/grace-acs-risk-mortality-calculator
HEART score	Score to assess risk of major adverse cardiac events (MACE) in ED patients with chest pain based on the following 5 predictors: History, ECG, Age, Risk factors and Troponin www.heartscore.nl
HF	heart failure

hs-cTn	high-sensitivity cardiac troponin
IHD	ischemic heart disease
kDa	kilo Dalton (unit of molecular weight)
LBBB	left bundle branch block
LoD	limit of detection
MACE	major adverse cardiac events
MI	myocardial infarction
MINOCA	myocardial infarction with non-obstructive coronary arteries
NPV	negative predictive value
NSTE-ACS	non-ST-segment elevation acute coronary syndromes
NSTEMI	non-ST-segment elevation myocardial infarction
NT-proBNP	N-terminal pro B-type natriuretic peptide
PCI	percutaneous coronary intervention
PE	pulmonary embolism
PPV	positive predictive value
RVD	right ventricular dysfunction
STEMI	ST-segment-elevation myocardial infarction
TIMI score	Thrombolysis in Myocardial Infarction; the TIMI risk score has been developed to assess mortality risk in patients with ACS www.mdcalc.com/timi-risk-score-ua-nstemi www.mdcalc.com/timi-risk-score-stemi
UA	unstable angina
URL	upper reference limit



bioMérieux

In vitro diagnostics serving public health

A major player in *in vitro* diagnostics for more than 50 years, bioMérieux has always been driven by a pioneering spirit and unrelenting commitment to improve public health worldwide.

Our diagnostic solutions bring high medical value to healthcare professionals, providing them with the most relevant and reliable information, as quickly as possible, to support treatment decisions and better patient care.

bioMérieux's mission entails a commitment to support medical education, by promoting access to diagnostic knowledge for as many people as possible. Focusing on the medical value of diagnostics, our collection of educational booklets aims to raise awareness of the essential role that diagnostic test results play in healthcare decisions.

Other educational booklets are available.
Consult your local bioMérieux representative,
or visit www.biomerieux.com/en/education/educational-booklets

The information in this booklet is for educational purposes only and is not intended to be exhaustive. It is not intended to be a substitute for professional medical advice. Always consult a medical director, physician, or other qualified health provider regarding processes and/or protocols for diagnosis and treatment of a medical condition. bioMérieux assumes no responsibility or liability for any diagnosis established or treatment prescribed by the physician.

bioMérieux S.A. • 69280 Marcy l'Étoile • France
Tel.: + 33 (0)4 78 87 20 00 • Fax: +33 (0)4 78 87 20 90
www.biomerieux.com